
ACCESSCUBE

Navigate with ease.

**GEBRAUCHSANWEISUNG
USER MANUAL**

Patented Product
developed by Medical Templates
www.medicaltemplates.ch

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CONTENTS

1 – Introduction	15
1.1 – Intended purpose	15
1.2 – Intended use	15
1.3 – Combination with other devices	16
1.4 – Overview of the device	16
2 – Important notes	18
2.1 – Meaning of the symbols	18
2.2 – Safety instructions	19
2.3 – Serious Incidents	19
2.4 – Patient safety Instructions	19
2.5 – Manual assembly	20
3 – Application	21
3.1 – Preparation	21
3.2 – Attaching the ACCESS CUBE	21
3.3 – Planning the puncture	22
3.4 – Placing the puncture needle	22
3.5 – Clearing the needle	23
3.6 – Performing the treatment	24
3.7 – Removing the ACCESS CUBE	24

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USER MANUAL

Thank you for purchasing our device. The ACCESS CUBE is characterized by its efficiency and ease of use. Please read the instruction manual carefully to ensure that you get the most out of the ACCESS CUBE. Please contact MEDICAL TEMPLATES or your local medical device adviser if you have any additional questions about this device or how to use it.

Manufacturer information

MEDICAL
TEMPLATES

MEDICAL TEMPLATES AG
Brunnenwiese 10
8132 Egg bei Zurich, Switzerland
info@medicaltemplates.ch
www.medicaltemplates.ch

1 – INTRODUCTION

This user manual forms an integral part of the device. It describes its intended use and must be kept near the device. Following the ACCESS CUBE instruction manual is essential for handling the device properly and making sure it works correctly. It ensures the safety of the patient and the user alike.

1.1 – Intended purpose

The ACCESS CUBE is used for guiding needles and instruments during Computed Tomography (CT)-guided punctures, such as infiltrations, diagnostic and therapeutic biopsies, or for placing access points for drainages or ablations. The ACCESS CUBE can be used anywhere on the body, provided that the ACCESS CUBE is resting on all four feet and is properly adhered to the patient. The device must be used in combination with a CT and associated software, by a physician within the radiology department. The ACCESS CUBE is a sterile, single-use device.

1.2 – Intended use

Medical indication

The ACCESS CUBE is used for guiding needles and instruments during CT-guided punctures, such as infiltrations, diagnostic and therapeutic biopsies, or for placing access points for drainages or ablations.

Patient profile

Adults and children from 3 years and up. It can be used anywhere on the body, provided that the ACCESS CUBE is resting on all four feet and is properly adhered to the patient.

User profile




The ACCESS CUBE can be unpacked and attached to the patient by a radiographer. A trained physician must perform all other functions.

1.3 – Combination with other devices

The following cannula and needle diameters are permitted for use with the ACCESS CUBE: G10 – G20.

1.4 – Overview of the device

A	Blister	1	Left half of the frame
B	ACCESS CUBE	2	Right half of the frame
C	Sealing film	3	Comb 1
		4	Comb 2
		5	Adhesive tab
		6	Peel-off film

2 – IMPORTANT NOTES

Important text passages that require special attention are marked with a warning triangle.

2.1 – Meaning of the Symbols

	Manufacturer
	Lot number
	Reference number
	Date of manufacture
	Use by date
	This device complies with the essential requirements set out in Directive 93/42/EEC concerning medical devices
	Consult instructions for use
	Sterilized using ethylene oxide
	Keep away from sunlight
	Keep dry
	Storage temperature: +10°C to +27°C
	Do not use if package is damaged
	Do not re-use
	Authorized representative in the European Community
	Importer in the European Community
	Medical device

2.2 – Safety instructions

The ACCESS CUBE is packed and shipped in assembled condition.

	The ACCESS CUBE is a sterile, single-use device. The ACCESS CUBE may only be handled when wearing sterile disposable gloves.
	Do not use the ACCESS CUBE if the packaging is damaged. Resterilization is not permitted, as the proper function of the Access Cube can no longer be guaranteed.
	The ACCESS CUBE may only be used in a sterile condition. Reuse jeopardizes the safety and health of the patient since the functionality of the adhesive tabs can no longer be guaranteed.

2.3 – Serious Incidents

Please report any event in connection with the ACCESS CUBE which presumably is a serious incident to Medical Templates AG (support@medicaltemplates.ch) and your local responsible authority. According to regulation (eu) 2017/745 on medical devices (MDR) a serious incident is any incident that directly or indirectly led, might have led or might lead to any of the following:

- (a) the death of a patient, user or other person
- (b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health
- (c) a serious public health threat

2.4 – Patient Safety Instructions

The ACCESS CUBE is affixed to the patient with adhesive tabs. The adhesive, Acrylate, is intended for medical use, but has a slight risk of causing minor skin irritation.

Never rely on the ACCESS CUBE alone for a puncture, as the target area can move after the initial scan due to breathing. Always perform a control scan to ensure correct positioning of the needle.


Always ensure the needle is inserted through the holes and corners indicated by the software. Using the wrong coordinates can lead to a missed target.

When inserting a needle or thin instrument through the ACCESS CUBE, there is a small risk that part of the ACCESS CUBE could be scraped by the needle and inserted into the patient. Place the needle carefully to avoid this.

2.5 – Manual assembly

The ACCESS CUBE should arrive fully assembled. If the combs shift during transportation, the following steps must be followed for manual assembly:

- 1 Hold comb 1 by the handle marked '1'.
- 2 With the other hand, hold comb 2 by the handle marked '2'.
- 3 From above, guide the tips of comb 1 at an angle into the corresponding spaces of comb 2.
- 4 Rest comb 1 on comb 2. The tips of the combs fit underneath the edge of the frame.
- 5 Place the joined combs into the top or bottom frame.

3 – APPLICATION

3.1 – Preparation

- 1 Identify the intervention site on the patient with the help of a topogram if needed and mark the location where the ACCESS CUBE is to be attached.
- 2 Prepare the location where the ACCESS CUBE is to be attached: Clean and disinfect the site, as well as the area surrounding the intervention site. Ideally, also place sterile covers around the ACCESS CUBE.
- 3 Carefully pull the sealing film off of the packaging and set it aside.

Do not dispose of the sealing film. It will be needed again at a later stage.

- 4 Wearing sterile gloves, remove the sterile ACCESS CUBE from its packaging.

3.2 – Attaching the ACCESS CUBE

- 1 Pull off the protective film covers off the adhesive tabs.
- 2 Position the ACCESS CUBE over the patient. The ACCESS CUBE can be mounted transversally or longitudinally.
- 3 Place the ACCESS CUBE on the patient's body in the desired location and press the ACCESS CUBE gently.
- 4 Press the adhesive tabs firmly against the skin to ensure the ACCESS CUBE is securely attached. Insufficient fixation of the feet can lead to the ACCESS CUBE shifting after the planning scan and an inaccurate puncture.

Once the adhesive tabs have been pressed on, repositioning is no longer possible.

3.3 – Planning the puncture

- 1 Perform a topogram to identify the scan area required for the puncture. Be sure it includes the ACCESS CUBE and puncture site.
- 2 Create a spiral CT image of the ACCESS CUBE.

The entire ACCESS CUBE must be captured. The ACCESS CUBE must not be covered or taped over at any point during the CT scan. Only then is identification of the ACCESS CUBE by planning software ensured.

- 3 Analyze the recorded images in compatible planning software to register the ACCESS CUBE and plan the puncture. Please refer to the planning software instruction manual for this purpose.
- 4 Note down the coordinates provided by the planning software on the form provided for this purpose, which can be found on the sealing film. Be sure to note which corner of the hole is indicated for the chosen route.

Ensure that the ACCESS CUBE's position has not changed between the planning scan and the puncture.

3.4 – Placing the puncture needle

- 1 Locate the coordinates specified by the software on the top frame of the ACCESS CUBE, with the help of the guidelines A – C and 1 – 3.

The markings next to the grids in the form of triangle, semicircle, circle and square are only used for identification by the software and are not used for puncture.

- 2 Position the puncture needle in the selected star-shaped hole, making sure the needle is also resting in the corner of the hole specified in the software.
- 3 Position the needle in the bottom frame using the same principle. Here, the guidelines contain the letters D – F and numbers 4 – 6.
- 4 Local anesthesia can now be administered at the injection site. To do this, place an anesthetic needle through any two holes in the ACCESS CUBE to reach the puncture site.

- 5 Now perform the planned puncture with the positioned needle. Ensure the needle is inserted through the holes and corners indicated by the software. We recommend inserting the needle 2 – 3 centimeters before performing a control scan. This allows small corrections of the needle angle, should this be necessary.

Do not apply pressure to the ACCESS CUBE while inserting the needle to ensure you do not change the position of the ACCESS CUBE. This is the only way to ensure that the needle position calculated by the planning software can be maintained.

While inserting the needle, the pressure on the needle must not be so high that it causes the needle to bend. Only apply pressure in the direction of the needle tip.

When using beveled needles, make sure that the needle follows a straight trajectory (e.g. rotate the beveled needle while piercing).

Finish the first puncture at least 2 cm before the proximal end of the needle meets the top combs. Then the combs can be removed without causing the needle to move, and the needle can be inserted deeper.

- 6 Perform a control scan to check the needle's position. Repeated control scans may be necessary to insert the needle safely.

3.5 – Clearing the needle

If necessary, the ACCESS CUBE can be completely removed so the needle can be moved freely or so other needles can be placed in the vicinity of the ACCESS CUBE. If doing so is not necessary, continue to »3.6 – Performing the treatment«.

To completely remove the ACCESS CUBE:

- 1 On the upper level, hold comb 2 steady by pressing the handle number '2'.
- 2 Now grasp comb 1 by the handle and tilt it upwards to release it from comb 2.
- 3 Remove comb 1 by moving it away from the cube. Be careful not to touch the needle.
- 4 Then remove comb 2, being careful not to touch the needle.
- 5 Repeat steps 1 – 4 for the bottom layer.
- 6 Carefully separate the right and left halves of the frame by first disconnecting the connections on the top and bottom levels one at a time. Start with the top level. Make sure to carefully detach the adhesive tabs from the patient's skin.

- 7 Remove the frame halves.
- 8 Now advance the needle until the target is reached. Additional control scans may be necessary.

3.6 – Performing the treatment

Perform the necessary treatment (e.g. injection, aspiration or ablation).

3.7 – Removing the ACCESS CUBE

- 1 Carefully remove the needle after treatment.
- 2 Carefully detach the entire ACCESS CUBE from the patient's skin.
- 3 Please dispose of the ACCESS CUBE and consumables you have used during the procedure according to local regulations and your facilities procedures.

Reuse of the ACCESS CUBE is not permitted.

The ACCESS CUBE could come into contact with blood or tissue during its use, especially when removing the needle from the patient. Therefore it should be regarded as contaminated with potentially infectious substances of human origin and handled with appropriate caution.

MEDICAL
TEMPLATES

MEDICAL TEMPLATES AG

Brunnenwiese 10
8132 Egg
Schweiz / Switzerland

info@medicaltemplates.ch
www.medicaltemplates.ch